CORMAY AMMONIA / ETHANOL CONTROLS

CORMAY AMMONIA / ETHANOL CONTROLS

- Nr kat. 5-163 3 x 3 ml
- AMMONIA / ETHANOL CONTROL 1 1 x 3 ml
 - AMMONIA / ETHANOL CONTROL 2 1 x 3 ml
 - AMMONIA / ETHANOL CONTROL 3 1 x 3 ml

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INTENDED USE

CORMAY AMMONIA / ETHANOL CONTROLS kit is intended for quality control of ethanol and ammonia determinations.

Controls are prepared from bovine serum albumin with preservatives and stabilizers added.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

STABILITY AND STORAGE

- Unopened controls are stable at 2-8°C until the expiration date printed on the package.
- Stability after opening: 20 days at 2-8°C tightly capped.

PROCEDURE

Controls are ready to use.

- Allow the controls to equilibrate to room temperature prior to use.
- Mix controls gently and accurately before use.
- Do not allow controls to stand for long time periods at room temperature. After use store vials tightly closed in refrigerator.

WARNINGS AND NOTES

- CORMAY AMMONIA / ETHANOL CONTROLS meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Ingredients:

CORMAY AMMONIA / ETHANOL CONTROLS contain reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H -isothiazol-3-one, mixture (3:1).

Warning:

H317 May cause an allergic skin reaction.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P302+P352 IF ON SKIN: Wash with plenty of soap and water.

P333+P313 If skin irritation or rash occurs: Get medical advice.

P363 Wash contaminated clothing before reuse.

- Do not use controls after expiry date
- Do not interchange the caps and plugs.
- Turbidity of controls may signify a contamination. Do not use turbid controls.
- Controls are not intended for use as a standard.

(EN) ASSIGNED VALUES

AMMONIA / ETHANOL CONTROL 1		Lot: 904-44
COMPONENT	MEAN	RANGE
Ammonia	69.8 µg/dl	49.4 – 90.3
	41.0 µmol/l	29.0 – 53.0
Ethyl alcohol	42.0 mg/dl	34.0 – 50.0
	9.11 mmol/l	7.38 – 10.9
AMMONIA / ETHANOL CONTROL 2		Lot: 904-44
COMPONENT	MEAN	RANGE
Ammonia	184 µg/dl	147 – 221
	108 µmol/l	86.0 – 130
Ethyl alcohol	103 mg/dl	82.0 – 124
	22.4 mmol/l	17.8 – 26.9
AMMONIA / ETHANOL CONTROL 3		Lot: 904-46
COMPONENT	MEAN	RANGE
Ammonia	465 µg/dl	371 – 559
	273 µmol/l	218 – 328
Ethyl alcohol	256 mg/dl	205 – 307
	55.6 mmol/l	44.5 – 66.6

Assigned values are specific for lot number of control material. It is recommended that each laboratory establish its own means and acceptable ranges.

WASTE MANAGEMENT

Please refer to local legal requirements

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