CORMAY HbA_{1c} DIRECT CALIBRATORS

CORMAY HbA_{1c} DIRECT CALIBRATORS (PL)

- Nr kat. 4-308 4 x 0,5 ml
- HbA_{1c} DIRECT CALIBRATOR 1 1 x 0,5 ml
 - HbA_{1c} DIRECT CALIBRATOR 2 1 x 0,5 ml
 - HbA_{1c} DIRECT CALIBRATOR 3 1 x 0,5 ml
 - HbA_{1c} DIRECT CALIBRATOR 4 1 x 0,5 ml

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INTENDED USE

CORMAY HbA_{1c} DIRECT CALIBRATORS kit is dedicated for preparation of the calibration curve for direct method HbA_{1c} determinations.

The calibrators must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

The determination of HbA_{1c} is most commonly performed for the evaluation of glycaemic control in diabetes mellitus. HbA_{1c} values provide an indication of glucose levels over the preceding 4-8 weeks.

RECONSTITUTION

Carefully remove the cap from vial and add exactly 0.5 ml of deionized water. Recap, gently mix for 10 minutes until all material has dissolved.

Before the assay HbA_{1c} DIRECT reconstituted calibrators should be treated with HEMOLYSING REAGENT according to HbA_{1c} DIRECT kit insert.

STABILITY AND STORAGE

1. Unopened calibrator is stable at 2-8°C until the expiration date printed on the vial. Protect from light and heat!
2. The constituent in reconstituted calibrator is stable for: 4 weeks at 2-8°C or 8 weeks at -20°C.

WARNINGS AND NOTES

- The calibrator is based on human blood that has been tested for HBsAg and antibodies to HIV and HCV and found to be non-reactive. However this material should be handled as though capable of transmitting infectious disease.
- Avoid contact with skin and mucous membranes.

ASSIGNED VALUES

Values in NGSP units [%]:

CORMAY HbA _{1c} DIRECT CALIBRATORS		Lot: 008-34			
Analyser	Calibrator				
	1	2	3	4	
Prestige 24i Biolis 24i Premium Biolis 15i Biolis 50i ACCENT-200 ACCENT-200 II GEN ACCENT-220S	5.9%	8.5%	12.5%	15.5%	
BS-400, BS-480 Olympus AU400/AU640 BS-800, BS-800M ACCENT-300 HITACHI 911/912	5.0%	7.9%	11.8%	15.3%	

TRACEABILITY

HbA_{1c} DIRECT calibrators are traceable to IFCC standards through the Tosoh calibrator.

WASTE MANAGEMENT

Please refer to local legal requirements.

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