CORMAY HbA_{1c} DIRECT CONTROLS

(PL)

CORMAY HbA_{1c} DIRECT CONTROLS

Nr kat. 4-328 4 x 0,5 ml

- HbA_{1c} DIRECT CONTROL N 2 x 0,5 ml
- HbA_{1c} DIRECT CONTROL P 2 x 0,5 ml

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CORMAY HbA_{1c} DIRECT CONTROLS

Cat. No 4-328 4 x 0,5 ml

- HbA_{1c} DIRECT CONTROL N 2 x 0.5 ml
- HbA_{1c} DIRECT CONTROL P 2 x 0.5 ml

INTENDED USE

CORMAY HbA_{1c} DIRECT CONTROLS (CONTROL N and CONTROL P) are designed for routine quality controls of haemoglobin A_{1c} concentration.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

The determination of HbA_{1c} is most commonly performed for the evaluation of glycaemic control in diabetes mellitus. HbA_{1c} values provide an indication of glucose levels over the preceding 4-8 weeks.

RECONSTITUTION

Open the vial carefully (CONTROL N or CONTROL P) and add exactly 0.5 ml of deionized water. Recap, gently mix for 10 minutes until all material has dissolved. The reconstituted controls may be dispensed and frozen at -20°C. Do not freeze and thaw more than once.

Before the assay HbA_{1c} DIRECT reconstituted controls should be treated with HEMOLYSING REAGENT according to HbA_{1c} DIRECT kit insert.

STABILITY AND STORAGE

1. Unopened control is stable at 2-8°C until the expiration date printed on the vial. Protect from light and heat!
2. Stability after reconstitution:
 - 1 month at 2-8°C
 - 2 months at -20°C.

WARNINGS AND NOTES

- The controls are based on human blood that has been tested for the HbsAg and antibodies to HIV and HCV and found to be non-reactive. However this material should be handled as though capable of transmitting infectious disease.
- Avoid contact with skin and mucous membranes.

ASSIGNED VALUES

Values in NGSP units [%]:

CONTROL	LOT NO *	TARGET [% HbA _{1c}]	RANGE [% HbA _{1c}]	ANALYSER
N	211-39	5,9	4,9 – 6,9	ACCENT-200 II GEN ACCENT-220S
P		12,1	10,6 – 13,6	
N	211-39	5,5	4,5 – 6,5	ACCENT-300 BS-400, BS-480 BS-800, BS-800M ACCENT S120 ACCENT M320 ACCENT M320 Biolis 24i Premium Biolis 30i
P		11,7	10,2 – 13,2	

* Revision: 05.09.2023

WASTE MANAGEMENT

Please refer to local legal requirements.

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