

CORMAY IMMUNO-CONTROL II (L)

Nr kat. / Cat. No 4-290

Seria / Lot: 902-54

Data ważn. / Exp.: 2019.12

PARAMETR / COMPONENT	METODA / METHOD	WARTOŚĆ ŚRODKOWA / ASSIGNED VALUE	ZAKRES / RANGE
MIOGLOBINA MYOGLOBIN	TIA-wzmacniane lateksem Latex-enhanced TIA	94 ng/ml	80 – 108 ng/ml
FERRYTYNA FERRITIN	TIA-wzmacniane lateksem Latex-enhanced TIA	105 ng/ml	89 – 121 ng/ml
CALKOWITE IgE TOTAL IgE	TIA-wzmacniane lateksem Latex-enhanced TIA	67 IU/ml	50 – 84 IU/ml

CORMAY IMMUNO-CONTROL II (L)

Nr kat. / Cat. No 4-290

Seria / Lot: 902-54

Data ważn. / Exp.: 2019.12

PARAMETR / COMPONENT	METODA / METHOD	WARTOŚĆ ŚRODKOWA / ASSIGNED VALUE	ZAKRES / RANGE
MIOGLOBINA MYOGLOBIN	TIA-wzmacniane lateksem Latex-enhanced TIA	94 ng/ml	80 – 108 ng/ml
FERRYTYNA FERRITIN	TIA-wzmacniane lateksem Latex-enhanced TIA	105 ng/ml	89 – 121 ng/ml
CALKOWITE IgE TOTAL IgE	TIA-wzmacniane lateksem Latex-enhanced TIA	67 IU/ml	50 – 84 IU/ml

CORMAY IMMUNO-CONTROL II (H)

Nr kat. / Cat. No 4-290

Seria / Lot: 902-55

Data ważn. / Exp.: 2019.12

PARAMETR / COMPONENT	METODA / METHOD	WARTOŚĆ ŚRODKOWA / ASSIGNED VALUE	ZAKRES / RANGE
MIOGLOBINA MYOGLOBIN	TIA-wzmacniane lateksem Latex-enhanced TIA	300 ng/ml	270 – 330 ng/ml
FERRYTYNA FERRITIN	TIA-wzmacniane lateksem Latex-enhanced TIA	396 ng/ml	337 – 455 ng/ml
CAŁKOWITE IgE TOTAL IgE	TIA-wzmacniane lateksem Latex-enhanced TIA	407 IU/ml	346 – 468 IU/ml

CORMAY IMMUNO-CONTROL II (H)

Nr kat. / Cat. No 4-290

Seria / Lot: 902-55

Data ważn. / Exp.: 2019.12

PARAMETR / COMPONENT	METODA / METHOD	WARTOŚĆ ŚRODKOWA / ASSIGNED VALUE	ZAKRES / RANGE
MIOGLOBINA MYOGLOBIN	TIA-wzmacniane lateksem Latex-enhanced TIA	300 ng/ml	270 – 330 ng/ml
FERRYTYNA FERRITIN	TIA-wzmacniane lateksem Latex-enhanced TIA	396 ng/ml	337 – 455 ng/ml
CAŁKOWITE IgE TOTAL IgE	TIA-wzmacniane lateksem Latex-enhanced TIA	407 IU/ml	346 – 468 IU/ml