

PROTEIN/ BIAŁKO	CALIBRATOR 1	CALIBRATOR 2	CALIBRATOR 3	CALIBRATOR 4	CALIBRATOR 5
ALBUMIN/ ALBUMINA	7.60 g/l	13.65 g/l	25.51 g/l	51.02 g/l	101.53 g/l
α -1 ANTITRYPSIN/ α -1 ANTYTRYPSYNA	0.20 g/l	0.42 g/l	0.89 g/l	1.80 g/l	3.67 g/l
α -1 GLYCOPROTEIN ACID/ KWAŚNA α 1-GLIKOPROTEINA	0.13 g/l	0.25 g/l	0.51 g/l	1.03 g/l	2.06 g/l
ANTITHROMBIN III/ ANTYTROMBINA III	0.038 g/l	0.076 g/l	0.175 g/l	0.362 g/l	0.724 g/l
COMPLEMENT C3/ SKŁADNIK C3 DOPEŁNIACZA	0.24 g/l	0.47 g/l	0.94 g/l	1.88 g/l	3.77 g/l
COMPLEMENT C4/ SKŁADNIK C4 DOPEŁNIACZA	0.043 g/l	0.085 g/l	0.170 g/l	0.340 g/l	0.680 g/l
CERULOPLASMIN*/ CERULOPLAZMINA*	0.095 g/l	0.165 g/l	0.303 g/l	0.596 g/l	1.220 g/l
HAPTOGLOBIN/ HAPTOGLOBINA	0.21 g/l	0.40 g/l	0.82 g/l	1.63 g/l	3.28 g/l
IgA	0.33 g/l	0.65 g/l	1.31 g/l	2.61 g/l	5.24 g/l
IgG	1.55 g/l	3.10 g/l	6.21 g/l	12.40 g/l	24.82 g/l
IgM	0.14 g/l	0.30 g/l	0.59 g/l	1.16 g/l	2.34 g/l
TRANSFERRIN/ TRANSFERYNA	0.40 g/l	0.79 g/l	1.58 g/l	3.16 g/l	6.32 g/l

*Values compatible with CERULOPLASMIN reagents since LOT number xxx-xx27

*Wartości kompatybilne z odczynnikami do oznaczania CERULOPLAZMINY serie od xxx-xx27

PROTEIN	CALIBRATOR 1	CALIBRATOR 2	CALIBRATOR 3	CALIBRATOR 4	CALIBRATOR 5
ALBUMIN/ ALBUMINA	7.60 g/l	13.65 g/l	25.51 g/l	51.02 g/l	101.53 g/l
α -1 ANTITRYPSIN/ α -1 ANTYTRYPSYNA	0.20 g/l	0.42 g/l	0.89 g/l	1.80 g/l	3.67 g/l
α -1 GLYCOPROTEIN ACID/ KWAŚNA α 1-GLIKOPROTEINA	0.13 g/l	0.25 g/l	0.51 g/l	1.03 g/l	2.06 g/l
ANTITHROMBIN III/ ANTYTROMBINA III	0.038 g/l	0.076 g/l	0.175 g/l	0.362 g/l	0.724 g/l
COMPLEMENT C3/ SKŁADNIK C3 DOPEŁNIACZA	0.24 g/l	0.47 g/l	0.94 g/l	1.88 g/l	3.77 g/l
COMPLEMENT C4/ SKŁADNIK C4 DOPEŁNIACZA	0.043 g/l	0.085 g/l	0.170 g/l	0.340 g/l	0.680 g/l
CERULOPLASMIN*/ CERULOPLAZMINA*	0.095 g/l	0.165 g/l	0.303 g/l	0.596 g/l	1.220 g/l
HAPTOGLOBIN/ HAPTOGLOBINA	0.21 g/l	0.40 g/l	0.82 g/l	1.63 g/l	3.28 g/l
IgA	0.33 g/l	0.65 g/l	1.31 g/l	2.61 g/l	5.24 g/l
IgG	1.55 g/l	3.10 g/l	6.21 g/l	12.40 g/l	24.82 g/l
IgM	0.14 g/l	0.30 g/l	0.59 g/l	1.16 g/l	2.34 g/l
TRANSFERRIN/ TRANSFERYNA	0.40 g/l	0.79 g/l	1.58 g/l	3.16 g/l	6.32 g/l

*Values compatible with CERULOPLASMIN reagents since LOT number xxx-xx27

*Wartości kompatybilne z odczynnikami do oznaczania CERULOPLAZMINY serie od xxx-xx27