

# CORMAY HbA1c DIRECT CONTROLS



KARTA WARTOŚCI WYZNACZONYCH / ASSAY SHEET

F07-QAD-01\_01 15.12.2023

Nr kat. / Cat. No	Seria / Lot	Data ważności / Expiration date
4-328	402-4343	2026-07

ANALIZATOR / ANALYSER	KONTROLA / CONTROL	WARTOŚĆ ŚRODKOWA / TARGET	ZAKRES / RANGE
ACCENT-200 II GEN ACCENT-220S	HbA1c DIRECT CONTROL N	<b>5.6</b> % <b>37.7</b> mmol/mol	<b>4.6 – 6.6</b> % <b>26.8 – 48.6</b> mmol/mol
	HbA1c DIRECT CONTROL P	<b>12.3</b> % <b>110.9</b> mmol/mol	<b>10.8 – 13.8</b> % <b>94.5 – 127.3</b> mmol/mol
ACCENT-300 BS-400, BS-480 BS-800, BS-800M ACCENT S120 ACCENT MC240 ACCENT M320 Biolis 24i Premium Biolis 30i	HbA1c DIRECT CONTROL N	<b>5.1</b> % <b>32.2</b> mmol/mol	<b>4.1 – 6.1</b> % <b>21.3 – 43.2</b> mmol/mol
	HbA1c DIRECT CONTROL P	<b>12.3</b> % <b>110.9</b> mmol/mol	<b>10.8 – 13.8</b> % <b>94.5 – 127.3</b> mmol/mol