



# Myt-CAL

HEMATOLOGY CALIBRATOR

**CALIBRATOR**

**LOT** PLUS0524



2024-06-05

Instrument	Parameter	Calibration Values	Tolerance Limits
<b>Orphée Mythic 18 II Generation</b>	WBC x 10 <sup>3</sup> /μL	9.0	± 0.2
	RBC x 10 <sup>6</sup> /μL	4.50	± 0.10
	HGB g/dL	13.5	± 0.2
	HCT%	38.5	± 1.0
	PLT x 10 <sup>3</sup> /μL	245	± 12
<b>Orphée Mythic 22</b>	WBC x 10 <sup>3</sup> /μL	8.9	± 0.2
	RBC x 10 <sup>6</sup> /μL	4.47	± 0.10
	HGB g/dL	12.2	± 0.2
	HCT%	36.9	± 1.0
	PLT x 10 <sup>3</sup> /μL	255	± 12
<b>Orphée Mythic 22 AL</b>	WBC x 10 <sup>3</sup> /μL	8.9	± 0.2
	RBC x 10 <sup>6</sup> /μL	4.44	± 0.10
	HGB g/dL	12.3	± 0.2
	HCT%	37.1	± 1.0
	PLT x 10 <sup>3</sup> /μL	252	± 12
<b>Orphée Mythic 60</b>	WBC x 10 <sup>3</sup> /μL	9.0	± 0.2
	RBC x 10 <sup>6</sup> /μL	4.44	± 0.10
	HGB g/dL	12.5	± 0.2
	MCV fL	80.0	± 2.0
	PLT x 10 <sup>3</sup> /μL	260	± 12
	MPV fL	9.6	± 0.5
<b>Orphée Mythic 70</b>	WBC x 10 <sup>3</sup> /μL	9.1	± 0.2
	RBC x 10 <sup>6</sup> /μL	4.48	± 0.10
	HGB g/dL	12.5	± 0.2
	MCV fL	81.5	± 2.0
	PLT x 10 <sup>3</sup> /μL	267	± 12
	MPV fL	10.5	± 0.5

All brands and products are trademarks or registered trademarks of their respective companies.

