

**HAEM 8 CONTROL L**  
**HAEM 8 CONTROL N**  
**HAEM 8 CONTROL H**



KARTA WARTOŚCI WYZNACZONYCH / ASSAY SHEET / АНАЛИТИЧЕСКА ЛИСТА

F07-QAD-01\_15 15.12.2023

Nr kat. / Cat. No / Номер кат.	Seria / Lot / Серия	Data ważności / Expiration date / Срок годности
<b>8-601</b>	<b>312-X0124</b>	<b>2024.04.05</b>
<b>8-602</b>		
<b>8-603</b>		

**Mindray BC 3000 PLUS**

(odczynniki / reagents / реагенты PZ CORMAY S.A.)

Parametr / Parameter / Параметр		HAEM 8 CONTROL					
		LOW		NORMAL		HIGH	
		Średnia / Mean / Средняя	Zakres ± / Range ± / Предел ±	Średnia / Mean / Средняя	Zakres ± / Range ± / Предел ±	Średnia / Mean / Средняя	Zakres ± / Range ± / Предел ±
WBC	10 <sup>9</sup> /l	2.0	± 0.5	10.0	± 1.0	19.7	± 2.5
RBC	10 <sup>12</sup> /l	2.38	± 0.18	4.75	± 0.24	5.86	± 0.30
Hb	g/dl	6.4	± 0.4	13.6	± 0.6	18.4	± 0.8
Hct	%	17.8	± 1.5	40.5	± 2.0	55.7	± 2.4
MCV	fl	74.9	± 5.0	85.3	± 5.0	95.2	± 5.0
MCH	pg	26.8	± 2.5	28.6	± 2.5	31.4	± 2.5
MCHC	g/dl	35.7	± 3.0	33.5	± 3.0	33.0	± 3.0
Plts	10 <sup>9</sup> /l	69	± 20	255	± 40	503	± 60
MPV	fl	8.7	± 3.0	8.4	± 3.0	8.5	± 3.0

**Mythic 18 Orphee**

(odczynniki / reagents / реагенты PZ CORMAY S.A.)

Parametr / Parameter / Параметр		HAEM 8 CONTROL					
		LOW		NORMAL		HIGH	
		Średnia / Mean / Средняя	Zakres ± / Range ± / Предел ±	Średnia / Mean / Средняя	Zakres ± / Range ± / Предел ±	Średnia / Mean / Средняя	Zakres ± / Range ± / Предел ±
WBC	10 <sup>9</sup> /l	2.2	± 0.5	10.0	± 1.0	20.2	± 2.5
RBC	10 <sup>12</sup> /l	2.28	± 0.18	4.50	± 0.24	5.61	± 0.30
Hb	g/dl	6.2	± 0.4	13.1	± 0.6	17.6	± 0.8
Hct	%	17.1	± 1.5	37.6	± 2.0	52.1	± 2.4
MCV	fl	75.1	± 5.0	83.6	± 5.0	92.8	± 5.0
MCH	pg	27.1	± 2.5	29.1	± 2.5	31.3	± 2.5
MCHC	g/dl	36.1	± 3.0	34.8	± 3.0	33.7	± 3.0
Plts	10 <sup>9</sup> /l	87	± 20	278	± 40	531	± 60
MPV	fl	8.8	± 3.0	8.5	± 3.0	8.6	± 3.0

**ABX: Micros (CBC-mode)**

(odczynniki / reagents / реагенты PZ CORMAY S.A.)

Parametr / Parameter / Параметр		HAEM 8 CONTROL					
		LOW		NORMAL		HIGH	
		Średnia / Mean / Средняя	Zakres ± / Range ± / Предел ±	Średnia / Mean / Средняя	Zakres ± / Range ± / Предел ±	Średnia / Mean / Средняя	Zakres ± / Range ± / Предел ±
WBC	10 <sup>9</sup> /l	2.1	± 0.5	10.3	± 1.0	21.0	± 2.5
RBC	10 <sup>12</sup> /l	2.27	± 0.18	4.60	± 0.24	5.76	± 0.30
Hb	g/dl	6.4	± 0.4	13.4	± 0.6	18.0	± 0.8
Hct	%	16.3	± 1.5	37.8	± 2.0	53.1	± 2.4
MCV	fl	72.0	± 5.0	82.0	± 5.0	92.0	± 5.0
MCH	pg	28.2	± 2.5	29.2	± 2.5	31.2	± 2.5
MCHC	g/dl	39.3	± 3.0	35.6	± 3.0	33.9	± 3.0
Plts	10 <sup>9</sup> /l	78	± 20	273	± 40	507	± 60
MPV	fl	10.1	± 3.0	9.4	± 3.0	9.3	± 3.0